- ESR Tubes without anticoagulant. Cat. No. 261270 (200 tubes) and 261170 (1000 tubes)
- ESR Tubes with anticoagulant. Cat. No. 261277 (200 tubes) and 261177 (1000 tubes)
- ESR Vacuum Tubes with anticoagulant. Cat. No. 265277 (200 tubes) and 265177 (1000 tubes)

**Intended Use** Venoplast® tubes are evacuated and non evacuated tubes designed for collection of whole blood by either venipuncture or manual blood transfer from primary blood tubes. The tubes are designed for quantitative testing of erythrocyte sedimentation rate (ESR) in venous whole blood samples using the Sediten® automated ESR reader and for use in a laboratory or point of care setting

**Summary and Explanation of the Test** Erythrocyte sedimentation rate (ESR) is defined as the rate of settlement of red blood cells in plasma in a vertical column in 1 hour. Citrated whole blood is added to a tube. The red blood cells rouleaux in the center of the tube and as they fall the plasma is displaced around them. Two methods were developed for the clinical laboratories: the Westergren and the Wintrobe methods. Both methods provide slightly different reference ranges and results. The Westergren method is more commonly used for determining ESR. The ESR performed on the Sediten instrument is an advancement over the manual one-hour method. Results of testing with the Sediten instrument correlate well with the traditional Westergren method, and testing is completed in 10 minutes.

**Principle of the Procedure** Venoplast vacuum tubes are designed to collect the sample by direct draw from the patient. Alternatively, EDTA whole blood may also be added to the Venoplast vacuum and non vacuum tubes manually by the user. The tubes are plastic and designed to minimize risk of exposure to blood samples. Once the mixed sample is placed in the tubes, the tubes are mixed and placed on the Sediten instrument. The instrument measures the rate of sedimentation of the erythrocytes in the sample, and calculates the patient result after 10 minutes.

**Reagents**

- Venoplast tubes are available with and without 0.32 ml of 3.8% Sodium Citrate and with and without vacuum
- For in vitro diagnostic use only.
- Venoplast tubes are stable through the labeled expiration date.

**Specimen Collection**

- Collect venous whole blood samples in Venoplast vacuum tubes by direct draw. If collecting the sample in another tube, manually transfer the sample to the Venoplast tube. The Sodium Citrate in the tubes should not be removed from the tube before collection by venipuncture or manual addition of sample. Add sample to the fill line marked on the tube. After adding the sample, mix the contents of the tube well using the Sedimix® ESR Tube Mixer.
- Venoplast tubes are designed to minimize potential exposure of the operator to blood fluids; however, leaks or spills resulting in exposure may occur. Use universal precautions for the prevention of the transmission of blood borne pathogens when handling blood samples.
- Follow state and local regulations and guidelines for disposal of biohazard waste.
- Interfering conditions inhibit or promote rouleau formation and result in falsely decreased or increased ESR results. Examples of these conditions are polycythemia, anemia, red blood cell abnormalities (acanthocytosis, poikilocytosis, anisocytosis, spherocytosis, sickle cell disease and microcytosis) and protein abnormalities such as hypofibrinogenemia, hypogammaglobulinemia and macroglobulinemia.

**Collected samples should be tested within 4 hours of collection, or up to 12 hours after collection if stored refrigerated. If stored refrigerated, collected samples should be allowed to come to room temperature and mixed thoroughly prior to analysis.**

An I-Button is included with each package of Venoplast tubes. The I-Button contains the inventory accounting system information of the Venoplast tubes for using with the Sediten ESR reader. The credit accounting system information must be loaded into the Sediten reader to enable the instrument to work.

**Instructions for use of the I-Button:**

1. Un-plug the Sediten Instrument.
2. While holding the I-Button against the I-Button probe, re-connect the Sediten to start the instrument.
3. On the readout of the Sediten, it can be observed that the number of credits being transferred to the unit (200 or 1000) depending on the package size of Venoplast tubes.
4. The credit accounting system will automatically adjust as the instrument is being used. When the inventory approaches zero, a new kit must be opened and the information from the I-Button loaded into the system.
5. When loading the credits from the I-Button, the instrument will recognize if there are credits remaining on the instrument and add the residual credits to the ones transferring from the I-Button.
6. The information from the I-Button may be transferred only once. Then the I-Button should be discarded.

**Procedure**

1. Mix filled sample with Sedimix ESR Tube Mixer before testing.
2. Insert the well mixed sample tubes into available channels on the Sediten. The instrument will begin timing automatically.
3. The testing cycle takes 10 minutes. After 10 minutes, the patient result will be displayed above the respective channel.
4. Remove tube from channel and discard. The channel is now ready for another test.

**Materials Provided** Venoplast Tubes I-Button Package Insert

**Materials Required but Not Provided**

- Sediten Automated ESR Reader, with Operator's Manual
- Sedimix ESR Tube Mixer (Recommended)

**Quality Control** Good laboratory practices recommend the use of appropriate controls. Polymedco offers Sed-Chek® 2 ESR controls at two levels, for use with the Sediten instrument. (Cat.no. ESR-2CT, ESR-2CTN, ESR-2CTA) Two levels of controls should be assayed at least once a day; however, each laboratory should follow their state and local requirements. Values obtained should fall within a specified range. If these values fall outside the range and repetition excludes error, contact Polymedco Technical Services at (800) 431-2123 extension 285 for assistance.

**Limitations of the Procedure**

1. Tubes are designed for single use only.
2. Allow samples that have been refrigerated to return to room temperature before testing.
3. Results may be affected by morphology of the red blood cells and protein abnormalities in the patient sample.
4. Ensure the Venoplast tubes are properly filled before testing. The Sediten will not perform testing on tubes filled over 2 mm above or under 6 mm below the fill line on the tube.
5. Ensure contents of tubes are adequately mixed before inserting them into the instrument. Venoplast tubes are designed for use solely with the Sediten ESR reader. Do not use with other automated ESR readers.

**Expected Values in mm/hr**

Age	Male	Female
to 50 years	0 – 15	0 – 20
over 50 years	0 – 20	0 – 30

**Performance Characteristics - Range of the Test**

Using the Venoplast tubes the Sediten can detect sedimentation rates from 0 – 140 mm/hr. Samples above 140 mm/hr will be reported as > 140.

**Precision**

Within Run Precision	Level 1	Level 2	Level 3
Mean (mm/ hr)	14.4	27.3	34.6
SD	1.51	2.60	3.46
% CV	10.5	9.6	10.0
n	8	8	8

Between Run Precision	Level 1	Level 2	Level 3
Mean (mm/ hr)	13.3	24.7	30.6
SD	1.16	2.48	2.22
% CV	8.7	10.0	7.2
n	20	20	20

**Correlation** This method (Y) was compared with another commercially available method (X) and the following linear regression was obtained:  $Y = 0.97X + 2.73$  and a correlation coefficient of  $r = 0.9810$ . 324 patient samples were analyzed spanning the range 1 mm/hr to 140 mm/hr.

**Bibliography**

- Henry, J.B., Clinical Diagnosis and Management by Laboratory Methods 20<sup>th</sup> ed. W.B.Saunders Co., Philadelphia (2001).
- Bridgen, M.L., Clinical Utility of the Erythrocyte Sedimentation Rate, Am Fam Physician, 1999 Oct 1; 60 (5): 1443-50.
- Bedell, S.E., Bust, B.T., Erythrocyte Sedimentation Rate, Am J Med 1985; 79: 1001-1009.
- Saadeh, C., Erythrocyte Sedimentation Rate: Old and New Clinical Application's, South Med J, 1968; 91 (3): 220-225.
- National Cancer Institute, U.S. National Institute of Health.

Manufactured by  
**LP ITALIANA SPA**  
 Via C. Reale, 15/4 - 20157 Milano (Italy)  
 Tel. +39 023933061 - Fax +39 0239313484  
 info@lpitaliana.com - www.lpitaliana.com

Distributed by  
**Polymedco, Inc.**  
 510 furnace Dock Road  
 Cortlandt Manor, NY 10567  
 Tel: 800-432-2121 or 914-739-5400 - Fax: 914-739-5890  
 info@polymedco.com-www.Polymedco.com

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- Henry, J.B., Clinical Diagnosis and Management by Laboratory Methods 20<sup>th</sup> ed. W.B. Saunders Co., Philadelphia (2001).
- Bridgen, M.L., Clinical Utility of the Erythrocyte Sedimentation Rate, Am Fam Physician, 1999 Oct 1; 60 (5): 1443-50.
- Bedell, S.E., Bust, B.T., Erythrocyte Sedimentation Rate. Am J Med 1985; 78: 1001-1009.
- Saadeh, C., Erythrocyte Sedimentation Rate: Old and New Clinical Applications, South Med J, 1998; 91 (3): 220-225.
- National Cancer Institute, U.S. National Institute of Health.

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