
Sediplast® Autozero Westergren ESR System	
Intended Use The Sediplast ESR system is for the determination of Erythrocyte Sedimentation Rate (ESR) in whole blood.	
Summary The ESR is a non-specific test which indicates the presence of inflammation in the body. The test is used as an initial screening tool and also as a follow-up test to monitor the effects of therapy and the progression or regression of disease.	
Principle A pipette (tube) is filled with 0.8 ml of anticoagulated blood (total blood column 200 mm) and allowed to stand for a period of one hour. After one hour, a reading is taken at the plasma meniscus. This reading represents the distance in mm from the zero mark on the tube and is equal to the Erythrocyte Sedimentation Rate.	
Components Provided 1) Fixed bore pipettes, with a self-zeroing cap, a reservoir for excess blood and a cap with a labyrinth path ending with a vented bore to allow the escape of the air. 2) Vials (which may or may not contain anticoagulant) and stopper caps with a predetermined fracture diaphragm at the bottom.	
General guidelines The Sediplast® system is intended for professional use only. All patient samples should be treated as potentially infectious and handled according to standard precautions. The Occupational Safety and Health Administration (OSHA) and Good Laboratory Practice recommend use of personal protection equipment (PPE) - such as protective eyewear, shields, gloves, lab coats and footwear - when working with or near potential hazardous materials. Preparation No component preparation is required. For those customers using vials without diluent, a 4:1 dilution of the whole blood is required. This is accomplished, for example, by mixing 0.8 ml of blood and 0.2 ml of 3.8% Sodium Citrate (or saline) in a clean disposable tube.	

INSTRUCTIONS FOR USE	
Assemble and check for integrity of supplies.	
1	Remove the pink stopper on the prefilled vial (0.2 ml of 3.8% Sodium Citrate is used as diluent). Using a transfer pipette, fill the vial to the BOTTOM of the indicated fill line, with 0.8 ml of blood to make the required 4:1 dilution ratio. The <i>Sediplast®</i> system is available with and without diluent.
2	Replace pink pierceable stopper and gently invert several times to obtain homogeneous mixing.
3	Place vial in its rack on a level counter surface free from vibration. The diaphragm of the pink stopper is calibrated to break under the light pressure made by inserting the pipette.
4	Carefully insert the pipette through the pierceable stopper and gently push downwards until the pipette comes into contact with the blood sample. Next, while holding the middle of the pipette, gently twist the pipette and push downwards until the pipette rests on the BOTTOM of the vial. DO NOT USE EXCESSIVE FORCE AND DO NOT HOLD OR BLOCK THE TOP OF THE PIPETTE
5	The pipette will autozero the blood and any REASONABLE excess will flow into the reservoir compartment.
6	To ensure proper results, it is essential that the pipette makes FIRM contact with the bottom of the vial.
7	Let sample stand for exactly one hour and then read the numerical results of the erythrocyte sedimentation in millimeters
8	Dispose of properly after use.
9	Remove gloves and wash hands in accordance with established procedures.
General guidelines	
10	<i>Testing Timetable</i> Blood at room temperature (15-25)°C: test should be set up within 4 hours Refrigerated blood: must first be brought to room temperature and thoroughly mixed before testing Blood kept at 4°C: test should be set up within 12 hours
11	<i>Normal Values</i> Male: (under 50) 0-15 mm/hr – (over 50) 0-20 mm/hr Female: (under 50) 0-20mm/hr – (over 50) 0-30mm/hr

FILLING INSTRUCTIONS

- 1) Remove the pink stopper on the prefilled vial (0.2 ml of 3.8% sodium citrate is used as diluent). Using a transfer pipette, fill the vial to the bottom of the indicated fill line with 0.8 ml of blood to make required 4:1 dilution.

Sediplast[®] system can be delivered also without diluent.

- 2) Replace pierceable stopper and gently invert several times to mix.
- 3) Place vial in its rack on a level surface. Carefully insert the pipette through the pierceable stopper until the pipette comes in contact with the bottom of the vial. The diaphragm of the pink stopper is calibrated to break under the light pressure made by inserting the pipette. The pipette will autozero the blood and any excess will flow into the reservoir compartment.
- 4) To ensure proper results, it is essential that the pipette makes **FIRM** contact with the bottom of the vial. Let sample stand for exactly one hour and then read the numerical results of the erythrocyte sedimentation in millimeters. Dispose of properly after use.

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